



**NATIONAL STANDARD OF THE
PEOPLE'S REPUBLIC OF CHINA**

GB 5009.120-2016

**National food safety standard - Determination of Sodium
Propionate and Calcium Propionate in Foods
(English Translation)**

Issued on: August 31, 2016

Implemented on: March 1, 2017

Issued by: National Health and Family Planning Commission of the PRC

Table of Contents

Foreword...	3
1 Scope...	4
2 Principle...	4
3 Reagents and materials...	4
4 Instruments and apparatuses...	5
5 Analysis steps...	5
6 Apparatus reference conditions...	6
7 Drawing of the standard curve...	6
8 Determination of sample solution...	7
9 Description of the analysis result...	7
10 Precision...	8
11 Others...	8
12 Principle...	8
13 Reagents and Materials...	8
14 Instruments and apparatuses...	9
15 Analysis steps...	9
16 Description of the analysis result...	10
17 Precision...	11
18 Others...	11
Appendix A Standard substance chromatogram...	12

1 Scope

This Standard specifies the determination methods of sodium propionate and calcium propionate in bean products, raw and moist flour products, bread, cakes, vinegar and soy sauce.

This Standard applies to the determination of sodium propionate and calcium propionate in bean products, raw and moist flour products, bread, cakes, vinegar and soy sauce.

Method 1 -- Liquid Chromatography

2 Principle

The propionate in the sample, through acidification, is converted into propionic acid; then, it is extracted by ultrasonic water bath or distilled by steam distilled; it is collected and adjusted in pH; it is determined by high performance liquid chromatography; and, the content of propionic acid is quantified by external standard method. The sodium propionate and calcium propionate in the sample are based on propionic acid; the contents of sodium propionate and calcium propionate can be calculated according to the respective parameters as needed.

3 Reagents and materials

Unless otherwise specified, all the reagents are analytical reagents, the water is grade-1 water which is specified by GB/T 6682.

3.1 Reagents

3.1.1

Phosphoric acid (H_3PO_4).

3.1.2

Diammonium hydrogen phosphate $[(\text{NH}_4)_2\text{HPO}_4]$.

3.1.3

Silicone oil.

3.2 Preparation of reagents

flour products, vinegar, soy sauce). after the sample is homogenized, accurately weigh 25 g (accurate to 0.01g) of it in a 500 mL distillation bottle; add 100 mL of water; then, use 50 mL of water to rinse the container; transfer it to a distillation flask; add 20 mL of 1 mol/L phosphoric acid solution, 2 ~ 3 drops of silicone oil for steam distillation; place the 250 mL volumetric flask in the ice bath as an absorption-solution device; when it is distilled to about 240 mL, take it out; leave it at room temperature for 30 min; use 1 mol/L phosphoric acid solution adjust the pH to about 3; add water to dilute to the scale; shake well; filter it through a 0.45 μm microporous membrane; use the liquid chromatography to measure later.

5.2.2 Direct digestion (applicable to bread, pastry). accurately weigh 5 g (accurate to

0.01 g) of sample into a 100 mL beaker; add 20 mL of water; add 0.5 mL of 1 mol/L phosphoric acid solution; mix well; after 10 min of ultrasonic extraction, use 1 mol/L phosphoric acid solution to adjust the pH to about 3; transfer the sample to a 50 mL volumetric flask; use water to fix-volume to the scale; shake well. Transfer all the samples to a 50 mL stoppered plastic centrifuge tube; centrifuge at not less than 4 000 r/min for 10 min; take the supernatant; after it is filtered through a 0.45 µm microporous membrane, place it aside for later liquid chromatographic determination.

6 Apparatus reference conditions

6.1

Chromatographic column. C18 column, 4.6 mm x 250 mm, 5 µm or equivalent column.

6.2 Mobile phase. 1.5 g/L diammonium hydrogen phosphate solution; use 1 mol/L

phosphoric acid solution to adjust the pH to 2.7~3.5 (prepare when needed); filter it through a 0.45 µm microporous membrane.

6.3 Flow velocity. 1.0 mL/min;

6.4 Column temperature. 25°C;

6.5 Injection volume. 20 µL;

6.6 Wavelength. 214 nm.

6.7 Reference conditions to clean the chromatographic column. after the experiment, use 10% methanol to wash for 1 h; then, use 100% methanol to wash for another 1 h.

7 Drawing of the standard curve

7.1

Distillation method. accurately absorb 0.5 mL, 1.0 mL, 2.5 mL, 5.0 mL, 7.5 mL,

10.0 mL and 12.5 mL of standard stock solutions in a 500 mL distillation flask; other operations are same as the sample preparation in 5.2.1; the final concentrations of its propionic acid standard solution are respectively 0.02 mg/mL, 0.04 mg/mL, 0.1 mg/mL, f -- dilution ratio.

The content of propionic acid that is measured in the sample multiply by the conversion coefficient of 1.296 7 equals the content of sodium propionate;

The content of propionic acid that is measured in the sample multiply by the conversion

coefficient of 1.256 9 equals the content of calcium propionate;

The calculation result shall keep three significant figures.

10 Precision

The absolute difference of two independent test results under repeatability cannot exceed 10% of the arithmetic mean value.

11 Others

When 25 g of sample is taken and the fix volume is 250 mL, the detection-limit of propionic acid is 0.03 g/kg, and the quantitation-limit is 0.10 g/kg.

Method 2 -- Gas chromatography

12 Principle

The propionate in the sample is converted to propionic acid by acidification; through steam distillation and collection, it directly enters the gas chromatography; use a hydrogen flame ionization detector to detect, so as to determine the retention time;

use the external standard method to quantify the content of propionic acid. The sodium propionate and calcium propionate in the sample are based on propionic acid; the contents of sodium propionate and calcium propionate can be calculated according to the respective parameters as needed.

13 Reagents and Materials

Note: Unless otherwise specified, all the reagents are analytical reagents, the water is grade-1 water which is specified by GB/T 6682.

13.1 Reagents

13.1.1

Phosphoric acid (H_3PO_4).

13.1.2

Formic acid (CH_2O_2).

absorption-solution device; when it is distilled to nearly 250 mL, take it out; leave it at room temperature for 30 min; add water to the scale. Mix well for gas chromatograph.