

ICS 27.120.01

CCS F 54



**NATIONAL STANDARD OF THE
PEOPLE'S REPUBLIC OF CHINA**

GB 13172-2009

**Molybdenum 99-Techetium 99m chromatographic generator
(fission)**

(English Translation)

Issued on: March 13, 2009

Implemented on: March 1, 2010

**Issued by: General Administration of Quality Supervision, Inspection and
Quarantine of the People's Republic of China**

Table of Contents

Foreword...	3
1 Scope...	5
2 Normative references...	5
3 Terms and definitions...	5
4 Product specifications...	6
5 Technical requirements...	7
6 Test methods...	8
7 Inspection rules...	10
8 Product marking, label and instruction manual...	11
9 Packaging, transportation and storage...	11

Foreword

All the technical contents of this standard are mandatory.

This standard replaces GB 13172-1991 "Molybdenum 99-Techne-99m chromatographic generator (fission)". As compared with GB 13172-1991, the main changes of this standard are as follows.

- a) MODIFY the standard name;
- b) ADD 5 normative references;
- c) CHANGE the test and determination in the original standard into inspection;
- d) DELETE clause 3.3 of the original standard; ADD such items as the elution, eluate, calibration time, and aluminum concentration, etc.;
- e) CHANGE the sodium chloride injection of the original standard into elution;
- f) DELETE the Chapter 4 from the original standard;
- g) ADD the new product specifications;
- h) DELETE the clause 6.1.1 and 6.1.4 from the original standard;
- i) ADD the technical requirements for radioactive contamination of the

generator surface;

j) ADJUST the clause 6.2.9 of the original standard into the technical requirements of the generator;

k) DELETE the clause 7.2 from the original standard;

l) ADD the chapter of inspection rules;

m) SPLIT the original chapter 8 into two chapters. Marking, sign and instruction manual, AND packaging, transportation, and storage.

This standard was proposed by China National Nuclear Corporation.

This standard shall be under the jurisdiction of the National Nuclear Energy Standardization Technical Committee.

The drafting organization of this standard. China Atomic Energy Research Institute.

The drafters of this standard. Niu Xuxiang, Chen Daming, Liu Guangyu.

Molybdenum 99-Techneium 99m

chromatographic generator (fission)

1 Scope

This standard specifies the technical requirements, inspection methods, inspection rules, marking, inspection certificate and instruction manual, packaging, transportation and storage of molybdenum 99-Techneium 99m chromatographic generator (fission).

This standard applies to the medical molybdenum 99-Techneium 99m chromatographic generator which uses the Al_2O_3 as the adsorbent; AND it is not applicable to the technetium 99m generator of other types.

2 Normative references

The provisions in following documents become the provisions of this Standard through reference in this Standard. For the dated references, the subsequent amendments (excluding corrections) or revisions do not apply to this Standard; however, parties who reach an agreement based on this Standard are

encouraged to study if the latest versions of these documents are applicable.

For undated references, the latest edition of the referenced document applies.

GB 11806 Regulations for the safe transport of radioactive material

GB/T 15239 Sampling procedures and tables for isolated lot inspection by attributes

EJ/T 843-1994 Radioactive nuclide activity measurement - Well-type ionization chamber method

EJ/T 845-1994 Radioactive drug solution - pH determination method

State Food and Drug Administration Note [2006] No.264 "Detailed rules for the radioactive drug instruction specification"

People's Republic of China Pharmacopoeia II (2005 version)

3 Terms and definitions

The following terms and definitions apply to this standard.

3.1

Molybdenum 99-technetium 99m chromatographic generator (fission)

^{89}Sr . $6 \times 10^{-5}\%$

^{90}Sr . $6 \times 10^{-6}\%$

alpha impurity nuclide. $1 \times 10^{-7}\%$

Other total beta, gamma nuclide. $1 \times 10^{-2}\%$

5.2.5 Radiochemical purity

$\text{Na}^{99\text{m}}\text{TcO}_4$ content shall not be less than 98%.

5.2.6 Activity

The activity of the $^{99\text{m}}\text{Tc}$ eluted from the generator shall be 90.0% ~ 110.0% of the indicated amount if calculated in accordance with the calibration time.

5.2.7 Bacterial endotoxin

It shall comply with the provisions of the People's Republic of China Pharmacopoeia II (2005 version).

5.2.8 Sterile

It shall comply with the provisions of the People's Republic of China Pharmacopoeia II (2005 version).

6.1.1 Eluting efficiency

In accordance with the calibration time and product specifications, CALCULATE the activity A1 of the 99mTc in the generator at the time of eluting; USE 10 mL of elution to elute the generator; DETERMINE the activity A2 of the 99mTc in the generator; USE the formula (1) to calculate the eluting efficiency η . Where.

η - Eluting efficiency, expressed as a percentage;

A1 - The activity of the 99mTc in the generator at the time of eluting, in becquerel (Bq);

A2 - The activity of the 99mTc in the eluate, in becquerel (Bq).

8.1 Marking

Generator outer packaging shall have the following markings.

- a) Radioactive marking;
- b) Radioactive transportation level marking;
- c) Markings of rainproof, handle with care, and face upwards.

8.2 Label

The label is divided into inner packaging label and outer packaging label. The inner packaging label is attached to the generator surface; AND the outer packaging label is attached to the outer packaging. The label shall have the following information.

- a) Product name;
- b) Radioactive marking;
- c) Approval document number;
- d) Specifications;