

ICS 65.080
CCS B 10

NY

**INDUSTRY STANDARD OF THE
PEOPLE'S REPUBLIC OF CHINA**

NY/T 2272-2012

**Soil amendment - Measurement of calcium, magnesium, and
silicon content**

(English Translation)

Issued on: December 24, 2012

Implemented on: January 1, 2013

Issued by: Ministry of Agriculture of the People's Republic of China

Table of Contents

Foreword...	3
1 Scope...	4
2 Normative references...	4
3 Measurement of calcium content...	4
4 Measurement of magnesium content...	10

Foreword

This standard was drafted in accordance with the provisions of GB/T 1.1.

This standard was proposed by AND shall be under the jurisdiction of the Ministry of Agriculture of the People's Republic of China.

The drafting organizations of this standard. National Fertilizer Quality Supervision and Inspection Center (Beijing).

1 Scope

This standard specifies the test methods for the soil conditioner - calcium, magnesium, AND silicon content measurement.

This standard applies to measurement of soil conditioner - calcium, magnesium, AND silicon content.

2 Normative references

The following documents are essential to the application of this document.

For the dated documents, only the versions with the dates indicated are applicable to this document; for the undated documents, only the latest version (including all the amendments) are applicable to this Standard.

HG/T 2843 Chemical fertilizer products - Standard volumetric, standard, reagent and indicator solutions for chemical analysis

NY/T 887 Density testing of liquid fertilizer

3.1.1 Principle

As for calcium of the sample solution, in slightly acidic medium, USE a certain amount of strontium salt as the release agent, to make it atomization in the lean combustion air - acetylene flame; the generated atomic vapor absorbs the characteristic wavelength 422.7 nm light which is emitted from the calcium hollow cathode lamp, AND absorbance values is directly proportional to the calcium ground state atomic concentration.

3.1.2 Reagents and materials

The preparation of the reagents, water and solution which are used in this standard shall, unless otherwise indicated of specifications AND preparation methods, be in accordance with the provisions of HG/T 2843.

It shall follow the provisions of 3.1.4.2.

3.2.4.3 Drawing of working curve

Respectively PIPETTE 0 mL, 0.50 mL, 1.00 mL, 2.00 mL, 4.00 mL, AND 5.00 mL of calcium standard solution (3.2.2.1) into six 100 mL volumetric flasks; USE water to make it reach to constant volume; MIX it uniformly. The calcium's mass concentration of this standard series solution is respectively 0 $\mu\text{g/mL}$, 5.0 $\mu\text{g/mL}$, 10.0 $\mu\text{g/mL}$, 20.0 $\mu\text{g/mL}$, 40.0 $\mu\text{g/mL}$, AND 50.0 $\mu\text{g/mL}$.

Before measurement, based on the nature of the element to be measured and the instrument performance, OPTIMIZE the measurement conditions of argon flow rate, observation height, RF generator power, AND integration time. Then USE the plasma emission spectrometer, at a wavelength of

317.933 nm, to measure the emission intensity of each standard solution.

USE the calcium's mass concentration of the standard series solution ($\mu\text{g/mL}$) as the abscissa AND the corresponding emission intensity as the ordinate; DRAW the working curve.

Note: It can, based on the different instrument sensitivities, ADJUST the mass

concentration of the standard series solution.

3.2.4.4 Measurement

USE the sample solution OR the solution which is diluted to certain ratio, under the same conditions of the standard series solution for measurement purposes, to measure the emission intensity of calcium; from the working curve, FIND the corresponding mass concentration of the calcium ($\mu\text{g/mL}$).

3.2.4.5 Blank test

Except for not adding of sample, the rest steps are same as those described in 3.2.4.4.

3.2.5 Analysis result presentation

It shall follow the provisions of 3.1.5.

3.2.6 Allowable difference

It shall follow the provisions of 3.1.6.

3.2.7 Conversion of mass concentration

It shall follow the provisions of 3.1.7.

4.1.4.1 Sample preparation

With respect to the solid samples, after multiple times of splitting, TAKE out about 100 g of sample; rapidly GRIND it to the grain which can all pass the

0.50 mm pore size sieve (if the sample is wet, it shall pass through the 1.00

mm pore size sieve); MIX it uniformly; PLACE it in a clean AND dry container; after the liquid sample is subject to multiple times of shaking, quickly TAKE out about 100 mL; PLACE it in a clean AND dry container.

4.1.4.2 Sample solution preparation

- Solid sample. WEIGH 0.2 g ~ 3 g of sample (accurate to 0.0001 g); PLACE into a 250 mL volumetric flask; ADD 150 mL of hydrochloric acid solution (4.1.2.5) which was pre-heated to 28 °C ~ 30 °C; tightly PLUG

the volumetric flask; SHAKE the volumetric flask to disperse the sample in the liquid; MAINTAIN the solution temperature at $28\text{ }^{\circ}\text{C} \sim 30\text{ }^{\circ}\text{C}$; USE the oscillator (4.1.3.2) the frequency of which is set to $(180 \pm 20)\text{ r/min}$ to oscillate it for 30 min; then TAKE the volumetric flask out; COOL it to room temperature; USE water to dilute it to the mark; MIX it uniformly; dry FILTER it; DISCARD the first few milliliters of filtrate; USE the rest filtrate for measurement.

- Liquid sample. WEIGH $0.2\text{ g} \sim 3\text{ g}$ sample (accurate to 0.0001 g); PLACE it into a 250 mL volumetric flask; USE water to make it reach to constant volume; MIX it uniformly; dry FILTER it; DISCARD the first few milliliters of filtrate; USE the rest filtrate for measurement.

4.1.4.3 Drawing of working curve

Respectively PIPETTE 0 mL, 1.00 mL, 2.00 mL, 4.00 mL, 8.00 mL, AND

10.00 mL of magnesium standard solution (4.1.2.4) into six 100 mL volumetric flasks; respectively ADD 4 mL of hydrochloric acid solution (4.1.2.1) AND 10 mL of strontium chloride solution (4.1.2.2); USE water to make it reach to constant volume; MIX it uniformly. The magnesium's mass concentration of this standard series solution is respectively $0\text{ }\mu\text{g/mL}$, $1.00\text{ }\mu\text{g/mL}$, $2.00\text{ }\mu\text{g/mL}$, $4.00\text{ }\mu\text{g/mL}$, $8.00\text{ }\mu\text{g/mL}$, AND $10.00\text{ }\mu\text{g/mL}$. Under the selected optimum working conditions, at the wavelength of 285.2 nm , USE the lean combustion air - acetylene flame; and USE the magnesium content

0 standard solution as the reference solution to set zero, in order to measure

the absorbance value of each standard solution. USE the magnesium's mass concentration of the standard series solution ($\mu\text{g/mL}$) as the abscissa AND the corresponding absorbance value as the ordinate; DRAW the working curve.